

Anti-Kdm4A Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - O75164

Item Number: CM0030979



Introduction

Rabbit polyclonal antibody against human KDM4A (JMJD2A), a histone demethylase that specifically removes methyl marks from H3K9me3 and H3K36me3. Validated for WB, IHC-P, IP, and ELISA with cross-reactivity to human and mouse. Suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	KDM4A Rabbit pAb
Remarks/Alias	JMJD2, JHDM3A, JMJD2A, TDRD14A
Species	Human
GeneID	9682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 585-815 of human KDM4A (NP_055478.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75164
Protein Weight	121kDa
Shipping	Ice bag