


Anti-Kdm4B/Jmjd2B Rabbit Polyclonal Antibody For Wb - O94953

Item Number: CM0029992



Introduction

High-quality rabbit polyclonal antibody against human KDM4B/JMJD2B, validated for WB and ELISA. Suitable for human and mouse samples. Ideal for epigenetic and neurological research.

[Learn More](#)

Parameter	Value
Product Name	KDM4B/JMJD2B Rabbit pAb
Remarks/Alias	MRD65; JMJD2B; TDRD14B; KDM4B/JMJD2B
Species	Human
GeneID (Human)	23030
GeneID	23030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-365 of human KDM4B/JMJD2B (NP_055830.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94953
Protein Weight	122kDa
Shipping	Ice bag