


Anti-Kdm4C Rabbit Polyclonal Antibody For Wb, Elisa - Q9H3R0

Item Number: CM0024717



Introduction

Rabbit polyclonal antibody targeting KDM4C (JMJD2C/GASC1) for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for epigenetic and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	KDM4C Rabbit pAb
Remarks/Alias	GASC1; JHDM3C; JMJD2C; TDRD14C; KDM4C
Species	Human
GeneID (Human)	23081
GeneID	23081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-830 of human KDM4C (NP_055876.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3R0
Protein Weight	120kDa
Shipping	Ice bag