

Anti-Kdm5A Polyclonal Antibody For Wb, Ip, Elisa, Chip - P29375

Item Number: CM0029209



Introduction

High-quality rabbit polyclonal antibody against human KDM5A (RBP2/RBBP2), validated for WB, IP, ELISA, and ChIP. Cross-reacts with human and rat. Targets lysine-specific demethylase 5A, a histone demethylase involved in chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	KDM5A Rabbit pAb
Remarks/Alias	RBP2; RBBP2; RBBP-2; KDM5A
Species	Human
GeneID	5927
Immunogen	Recombinant protein of human KDM5A
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA, ChIP
Cross Reactivity	Human, Rat
SWISS	P29375
Protein Weight	192kDa
Shipping	Ice bag