


Anti-Kdm5B Rabbit Polyclonal Antibody [Ko Validated] For Wb, If/Icc, Elisa - Q9Ugl1

Item Number: CM0029608



Introduction

Rabbit polyclonal anti-KDM5B antibody, KO-validated, for WB, IF/ICC, ELISA. Recognizes human and mouse KDM5B (Lysine-specific demethylase 5B). Ideal for chromatin modification and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] KDM5B Rabbit pAb
Remarks/Alias	CT31; PLU1; PUT1; MRT65; PLU-1; JARID1B; PPP1R98; RBP2-H1; RBBP2H1A; 5B
Species	Human
GeneID (Human)	10765
GeneID	10765
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 784-883 of human KDM5B (NP_006609.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UGL1
Protein Weight	176kDa
Shipping	Ice bag