


Anti-Kdm6A (Utx) Rabbit Mab [Kd Validated] For Wb, If-P, Ihc-P, Elisa - O15550

Item Number: CM0005411



Introduction

High-quality anti-KDM6A rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Recognizes endogenous human, mouse, rat KDM6A. Provides reliable detection of lysine-specific demethylase 6A (UTX) for epigenetic and developmental research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] KDM6A Rabbit mAb
Remarks/Alias	UTX; KABUK2; bA386N14.2
Species	Human
GeneID (Human)	7403
GeneID	7403
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 611-710 of human KDM6A (NP_066963.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15550
Protein Weight	154kDa
Shipping	Ice bag