


Anti-Kdm7A Polyclonal Antibody For Ihc-P And Elisa - Q6Zmt4

Item Number: CM0019132



Introduction

Rabbit polyclonal KDM7A antibody validated for IHC-P and ELISA. Reacts with human, mouse, and rat. 106kDa. Suitable for studying histone demethylation and chromatin regulation in neurodevelopment.

[Learn More](#)

Parameter	Value
Product Name	KDM7A Rabbit pAb
Remarks/Alias	JHDM1D; KDM7A
Species	Human
GeneID (Human)	80853
GeneID	80853
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 840-941 of human KDM7A (NP_085150.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZMT4
Protein Weight	106kDa
Shipping	Ice bag