

Anti-Kdm7A Rabbit Polyclonal Antibody For Wb And Elisa - Q3Uwm4

Item Number: CM0029883

Introduction

Rabbit polyclonal antibody targeting KDM7A, a histone demethylase crucial for brain development. Validated for WB and ELISA, with cross-reactivity in human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	KDM7A Rabbit pAb
Remarks/Alias	Jhdm1d; mKIAA1718; A630082K20Rik; KDM7A
Species	Mouse
GeneID (Human)	80853
GeneID (Mouse)	338523
Immunogen	Synthetic peptide or Recombinant fusion protein containing amino acids 417-735 of mouse KDM7A (NP_001028602.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3UWM4
Protein Weight	92kDa/106kDa
Shipping	Ice bag