

Anti-Kdm7A Rabbit Polyclonal Antibody For Wb, Elisa - Q3Uwm4

Item Number: CM0013136

Introduction

KDM7A Rabbit Polyclonal Antibody for Western blot (WB) and ELISA, cross-reactive with human. Lysine-specific demethylase 7A regulates histone methylation, essential for brain development and chromatin dynamics.

[Learn More](#)

Parameter	Value
Product Name	KDM7A Rabbit pAb
Remarks/Alias	Jhdm1d; mKIAA1718; A630082K20Rik; KDM7A
Species	Mouse
GeneID (Human)	80853
GeneID (Mouse)	338523
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of mouse KDM7A (NP_001028602.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q3UWM4
Protein Weight	106kDa
Shipping	Ice bag