


Anti-Kdm8 (Jmjd5) Rabbit Polyclonal Antibody For Wb, Elisa - Q8N371

Item Number: CM0016821



Introduction

KDM8 (JMJD5) rabbit polyclonal antibody for WB, ELISA. Cross-reacts with human and mouse. Ideal for chromatin biology, cell cycle, and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	KDM8 Rabbit pAb
Remarks/Alias	JMJD5; KDM8
Species	Human
Gene ID	79831
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 147-416 of human KDM8 (NP_079049.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N371
Protein Weight	47kDa
Shipping	Ice bag