

Anti-Keap1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14145

Item Number: CM0023531



Introduction

High-quality rabbit polyclonal antibody against KEAP1, validated for WB, IF/ICC, ELISA in human, mouse, and rat. Ideal for oxidative stress and NRF2 pathway research. SWISS: Q14145.

[Learn More](#)

Specification	Detail
Product Name	KEAP1 Rabbit pAb
Remarks/Alias	INrf2; KLHL19; KEAP1
Species	Human
Gene ID (Human)	9817
Gene ID	9817
Immunogen	Synthetic Peptide. A synthetic peptide corresponding to a sequence within amino acids 301-400 of human KEAP1 (NP_036421.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14145
Protein Weight	70kDa
Shipping	Ice bag