

Anti-Keap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14145

Item Number: CM0019198



Introduction

This rabbit polyclonal antibody targets human KEAP1 (Kelch-like ECH-associated protein 1), the oxidative stress sensor. Suitable for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	KEAP1 Rabbit pAb
Remarks/Alias	INrf2; KLHL19; KEAP1
Species	Human
GeneID	9817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-510 of human KEAP1 (NP_036421.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14145
Protein Weight	70kDa
Shipping	Ice bag