

Anti-Khdrbs1/Sam68 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q07666

Item Number: CM0029864



Introduction

Rabbit polyclonal antibody targeting human KHDRBS1/Sam68, suitable for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat. Key roles in signal transduction, alternative splicing, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	KHDRBS1/Sam68 Rabbit pAb
Remarks/Alias	p62; p68; Sam68; KHDRBS1/Sam68
Species	Human
Gene ID (Human)	10657
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07666
Protein Weight	48kDa
Shipping	Ice bag