

Anti-Khdrbs3 Polyclonal Antibody For Wb, If/Icc, Elisa - 075525

Item Number: CM0029446



Introduction

High-quality rabbit polyclonal antibody targeting human KHDRBS3, validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse. Ideal for research on alternative splicing regulation and RNA-binding protein functions.

[Learn More](#)

Parameter	Value
Product Name	KHDRBS3 Rabbit pAb
Remarks/Alias	Etle; SALP; SLM2; SLM-2; TSTAR; T-STAR; etoile; KHDRBS3
Species	Human
Gene ID	10656
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 177-346 of human KHDRBS3 (NP_006549.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75525
Protein Weight	39kDa
Shipping	Ice bag