

Anti-Khk Rabbit Polyclonal Antibody For Wb, Elisa - P50053

Item Number: CM0026073



Introduction

KHK Rabbit Polyclonal Antibody for WB, ELISA, reacts with human, mouse, and rat. Target is ketohexokinase (hepatic fructokinase) involved in fructose metabolism. Ideal for carbohydrate metabolism and kinase research.

[Learn More](#)

Parameter	Value
Product Name	KHK Rabbit pAb
Remarks/Alias	FRUCTU; KHK
Species	Human
GeneID (Human)	3795
GeneID	3795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-298 of human KHK (NP_000212.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50053
Protein Weight	33kDa
Shipping	Ice bag