


Anti-Ki67 Rabbit Polyclonal Antibody For If/Icc, Elisa - P46013

Item Number: CM0017811



Introduction

Rabbit polyclonal antibody targeting human Ki67 (P46013), a proliferation marker. Validated for IF/ICC and ELISA. Raised against a synthetic peptide (aa 700-800). Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Ki67 Rabbit pAb
Remarks/Alias	KiA, MIB-, MIB-1, PPP1R105, Ki67
Species	Human
GeneID (Human)	4288
GeneID	4288
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-800 of human Ki67 (NP_002408.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P46013
Protein Weight	359kDa
Shipping	Ice bag