

Anti-Ki67 Rabbit Polymab® (Ko Validated) For Wb, If/Icc, If-P, Ihc-P, Elisa - E9Pvx6

Item Number: CM0005030

Introduction

Ki67 Rabbit PolymAb® (KO validated), an essential tool for cell proliferation studies. Reacts with human, mouse, rat Ki67. Validated for WB, IF/ICC, IF-P, IHC-P, ELISA. ~350kDa protein. - E9PVX6

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Ki67 Rabbit PolymAb®
Remarks/Alias	KIA; MIB-; MIB-1; PPP1R105
Species	Human
Gene ID (Human)	17345
Gene ID	17345
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	E9PVX6
Protein Weight	350kDa/359kDa
Shipping	Ice bag