

Anti-Kiaa0196 (Washc5) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12768

Item Number: CM0024408

Introduction

Rabbit polyclonal antibody targeting human KIAA0196 (strumpellin, WASHC5), suitable for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for studying WASH complex-mediated endosomal actin dynamics and protein trafficking in IVD research.

[Learn More](#)

Parameter	Value
Product Name	KIAA0196 Rabbit pAb
Remarks/Alias	RTSC; SPG8; RTSC1; KIAA0196
Species	Human
GeneID	9897 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-530 of human KIAA0196 (NP_055661.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q12768
Protein Weight	134kDa
Shipping	Ice bag