

Anti-Kiaa1456 Polyclonal Antibody For Wb - Q9P272

Item Number: CM0024142

Introduction

Rabbit polyclonal antibody targeting KIAA1456/TRMT9B, a probable tRNA methyltransferase with tumor suppressor function. Validated for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for cancer and epitranscriptomics research.

[Learn More](#)

Parameter	Value
Product Name	KIAA1456 Rabbit pAb
Remarks/Alias	KIAA1456; C8orf79; TRM9L
Species	Human
GeneID	57604
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human KIAA1456 (NP_065895.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P272
Protein Weight	14 kDa
Shipping	Ice bag