

Anti-Kiaa1524 (Cip2A) Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tcg1

Item Number: CM0015031



Introduction

Rabbit polyclonal antibody targeting human KIAA1524 (CIP2A), validated for WB and ELISA. Suitable for PP2A inhibition, MYC stabilization studies, and genome instability research in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	KIAA1524 Rabbit pAb
Remarks/Alias	p90; NOCIVA; KIAA1524
Species	Human
GeneID (Human)	57650
GeneID	57650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 681-905 of human KIAA1524 (NP_065941.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TCG1
Protein Weight	102kDa
Shipping	Ice bag