

Anti-Kif13A Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H1H9

Item Number: CM0013709

Introduction

Rabbit polyclonal antibody against human KIF13A, validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Suitable for studying kinesin-mediated intracellular transport, melanosome biogenesis, and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	KIF13A Rabbit pAb
Remarks/Alias	RBKIN; bA500C11.2; KIF13A
Species	Human
GeneID (Human)	63971
GeneID	63971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1300-1600 of human KIF13A (NP_001099036.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H1H9
Protein Weight	202kDa
Shipping	Ice bag