

Anti-Kif14 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15058

Item Number: CM0015068

Introduction

Anti-KIF14 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC and ELISA. Detects human, mouse, rat KIF14. Suitable for cell division, cytokinesis, and apoptosis research. Ideal for studies of microtubule motor function.

[Learn More](#)

Parameter	Value
Product Name	KIF14 Rabbit pAb
Remarks/Alias	MKS12, MCPH20, KIF14
Species	Human
Gene ID	9928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1399-1648 of human KIF14 (NP_055690.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15058
Protein Weight	186kDa
Shipping	Ice bag