

Anti-Kif17 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9P2E2

Item Number: CM0019129

Introduction

High-quality rabbit polyclonal antibody against human KIF17, validated for WB, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat KIF17. Ideal for studying ciliary trafficking and neuronal transport.

[Learn More](#)

Parameter	Value
Product Name	KIF17 Rabbit pAb
Remarks/Alias	KIF3X; KLP-2; OSM-3; KIF17B; KIF17
Species	Human
GeneID (Human)	57576
GeneID	57576
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 360-660 of human KIF17 (NP_001116291).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2E2
Protein Weight	115kDa
Shipping	Ice bag