

Anti-Kif19 Rabbit Polyclonal Antibody For Wb, Elisa - Q2Tac6

Item Number: CM0018088

Introduction

KIF19 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat KIF19. Ideal for cilia and microtubule dynamics research.

[Learn More](#)

Parameter	Value
Product Name	KIF19 Rabbit pAb
Remarks/Alias	KIF19A; KIF19
Species	Human
GeneID (Human)	124602
GeneID	124602
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q2TAC6
Protein Weight	111kDa
Shipping	Ice bag