

Anti-Kif1C Polyclonal Antibody For Wb, Elisa - O43896

Item Number: CM0023759



Introduction

KIF1C Rabbit pAb is a polyclonal antibody validated for WB and ELISA applications, detecting human, mouse, and rat KIF1C. Ideal for studies of Golgi-to-ER retrograde transport and microtubule motor function.

[Learn More](#)

Parameter	Value
Product Name	KIF1C Rabbit pAb
Remarks/Alias	SAX2; LTXS1; SATX2; SPAX2; SPG58; KIF1C
Species	Human
GeneID (Human)	10749
GeneID	10749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 810-960 of human KIF1C (NP_006603.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43896
Protein Weight	123kDa
Shipping	Ice bag