

Anti-Kif20A Polyclonal Antibody For Wb, Elisa - O95235

Item Number: CM0015901



Introduction

KIF20A rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Detects a 100 kDa protein involved in mitosis and Golgi dynamics. Ideal for cell biology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	KIF20A Rabbit pAb
Remarks/Alias	RCM6; MKLP2; RAB6KIFL; KIF20A
Species	Human
GeneID (Human)	10112
GeneID	10112
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 591-890 of human KIF20A (NP_005724.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95235
Protein Weight	100kDa
Shipping	Ice bag