

Anti-Kif20B Rabbit Polyclonal Antibody For Wb, Elisa - Q96Q89

Item Number: CM0015843

Introduction

Rabbit polyclonal antibody targeting human KIF20B, specifically its C-terminal region. Validated for Western blot and ELISA applications; cross-reacts with mouse and rat. Suitable for studying cytokinesis, neuronal polarization, and bladder cancer.

[Learn More](#)

Parameter	Value
Product Name	KIF20B Rabbit pAb
Remarks/Alias	CT90, MPP1, KRMP1, MPP-1, MPHOSPH1, KIF20B
Species	Human
Gene ID (Human)	9585
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1200-1500 of human KIF20B (NP_057279.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Q89
Protein Weight	211kDa
Shipping	Ice bag