

Anti-Kif23 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q02241

Item Number: CM0009966



Introduction

Rabbit monoclonal antibody against human KIF23 (MKLP1), validated for Western blotting (WB), immunoprecipitation (IP), and ELISA. Recognizes human and rat KIF23, with a predicted molecular weight of ~110 kDa. Suitable for cytokinesis and cell division studies.

[Learn More](#)

Parameter	Value
Product Name	KIF23 Rabbit mAb
Remarks/Alias	CDA3; CHO1; CDAN3; KNSL5; MKLP1; CDAIII; CDAN3A; MKLP-1; KIF23
Species	Human
GeneID (Human)	9493
GeneID	9493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 801-960 of human KIF23 (Q02241)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q02241
Protein Weight	110kDa
Shipping	Ice bag