

Anti-Kif2A Polyclonal Antibody For Wb And Elisa - 000139

Item Number: CM0019989



Introduction

KIF2A rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Detects human, mouse, rat KIF2A (~80kDa). Suitable for mitosis and neuronal development research.

[Learn More](#)

Parameter	Value
Product Name	KIF2A Rabbit pAb
Remarks/Alias	HK2; KIF2; CDCBM3; KIF2A
Species	Human
Gene ID (Human)	3796
Gene ID	3796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-220 of human KIF2A (NP_004511.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00139
Protein Weight	80kDa
Shipping	Ice bag