


Anti-Kif2C/Mcak Rabbit Monoclonal Antibody For Wb, Elisa - Q99661

Item Number: CM0005145



Introduction

Rabbit monoclonal antibody targeting human KIF2C/MCAK (kinesin-like protein KIF2C, mitotic centromere-associated kinesin). Suitable for WB and ELISA. Research use in mitosis, chromosome segregation, and kinetochore-microtubule dynamics.

[Learn More](#)

Parameter	Value
Product Name	KIF2C/MCAK Rabbit mAb
Remarks/Alias	MCAK; CT139; KNSL6; KIF2C/MCAK
Species	Human
GeneID (Human)	11004
GeneID	11004
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIF2C/KIF2C/MCAK (Q99661).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99661
Protein Weight	81kDa
Shipping	Ice bag