

Anti-Kif3A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Y496

Item Number: CM0029767

Introduction

KIF3A Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, IP, ELISA applications. Detects KIF3A, a microtubule motor protein involved in primary cilium formation and intracellular transport. Reacts with Human, Mouse, Rat.

[Learn More](#)

Parameter	Value
Product Name	KIF3A Rabbit pAb
Remarks/Alias	FLA10; KLP-20; KIF3A
Species	Human
Gene ID (Human)	11127
Gene ID	11127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-699 of human KIF3A (NP_008985.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y496
Protein Weight	80kDa
Shipping	Ice bag