

Anti-Kif3A Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y496

Item Number: CM0012500

Introduction

Rabbit monoclonal antibody against human KIF3A, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying ciliary transport, centriole function, and hedgehog signaling.

[Learn More](#)

Parameter	Value
Product Name	KIF3A Rabbit mAb
Remarks/Alias	FLA10; KLP-20; KIF3A
Species	Human
GeneID (Human)	11127
GeneID	11127
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 600-699 of human KIF3A (Q9Y496).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y496
Protein Weight	80kDa
Shipping	Ice bag