


Anti-Kif4A Monoclonal Antibody For Wb - O95239

Item Number: CM0010712



Introduction

KIF4A Rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat KIF4A. Ideal for studying chromosome segregation and cytokinesis in mitosis.

[Learn More](#)

Parameter	Value
Product Name	KIF4A Rabbit mAb
Remarks/Alias	KIF4; KIF4G1; MRX100; XLID100; KIF4A
Species	Human
GeneID (human)	24137
GeneID	24137
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIF4A (O95239).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95239
Protein Weight	140kDa
Shipping	Ice bag