


Anti-Kif4A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95239

Item Number: CM0013490



Introduction

Rabbit polyclonal antibody targeting human KIF4A, a chromosome-associated kinesin essential for mitosis and cytokinesis. Validated for Western blot, IHC-P, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	KIF4A Rabbit pAb
Remarks/Alias	KIF4; KIF4G1; MRX100; XLID100; KIF4A
Species	Human
GeneID (Human)	24137
GeneID	24137
Immunogen	Recombinant protein corresponding to amino acids 870-1080 of human KIF4A (NP_036442.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95239
Protein Weight	140kDa
Shipping	Ice bag