


Anti-Kif5A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q12840

Item Number: CM0028962



Introduction

Anti-KIF5A Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA against human, mouse, rat. Targets neuron-specific kinesin heavy chain, applicable in neurobiology and axonal transport studies.

[Learn More](#)

Parameter	Value
Product Name	KIF5A Rabbit pAb
Remarks/Alias	NKHC; ALS25; MYO50; NEIMY; SPG10; D12S1889; KIF5A
Species	Human
GeneID (human)	3798
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 933-1032 of human KIF5A (NP_004975.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12840
Protein Weight	117kDa
Shipping	Ice bag