

Anti-Kif5B Monoclonal Antibody For Wb, If-P, Elisa - P33176

Item Number: CM0003445



Introduction

KIF5B Rabbit mAb validated for WB, IF-P, ELISA. Detects endogenous KIF5B in human, mouse, rat. Ideal for research on kinesin-1 motor function, mitochondrial/lysosomal transport, and axon elongation.

[Learn More](#)

Parameter	Value
Product Name	KIF5B Rabbit mAb
Remarks/Alias	KNS; KINH; KNS1; UKHC; HEL-S-61; KIF5B
Species	Human
GeneID (Human)	3799
GeneID	3799
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIF5B (NP_004512.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33176
Protein Weight	110kDa
Shipping	Ice bag