

Anti-Kif5B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P33176

Item Number: CM0017874



Introduction

KIF5B Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, and rat samples. Target kinesin-1 heavy chain for motor protein and intracellular transport research.

[Learn More](#)

Parameter	Value
Product Name	KIF5B Rabbit pAb
Remarks/Alias	KNS; KINH; KNS1; UKHC; HEL-S-61; KIF5B
Species	Human
GeneID (Human)	3799
GeneID	3799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 650-770 of human KIF5B (NP_004512.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33176
Protein Weight	110kDa
Shipping	Ice bag