

Anti-Kif7 Rabbit Polyclonal Antibody For Wb, Elisa - Q2M1P5

Item Number: CM0021489



Introduction

KIF7 rabbit polyclonal antibody for WB and ELISA applications. Raised against recombinant human KIF7 (aa1044-1343). Detects 151kDa protein. Cross-reacts with human. Ideal for hedgehog signaling studies.

[Learn More](#)

Parameter	Description
Product Name	KIF7 Rabbit pAb
Remarks/Alias	ACLS, AGBK, HLS2, MMEDF, JBTS12, UNQ340, KIF7
Species	Human
GeneID (human)	374654
GeneID	374654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1044-1343 of human KIF7 (NP_940927.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q2M1P5
Protein Weight	151kDa
Shipping	Ice bag