

Anti-Kifap3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q92845

Item Number: CM0010902

Introduction

KIFAP3 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets kinesin-associated protein 3, essential for chromosome movement and spindle tethering. Ideal for cell division and motor protein research.

[Learn More](#)

Parameter	Value
Product Name	KIFAP3 Rabbit mAb
Remarks/Alias	FLA3; KAP3; SMAP; KAP-1; KAP-3; Smg-GDS; dj190I16.1; KIFAP3
Species	Human
GeneID (Human)	22920
GeneID	22920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 675-792 of human KIFAP3 (Q92845)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92845
Protein Weight	91kDa
Shipping	Ice bag