

## Anti-Kifap3 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q92845

Item Number: CM0027228

### Introduction

Anti-KIFAP3 rabbit polyclonal antibody validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat KIFAP3 involved in chromosome movement. Ideal for cell cycle and mitosis research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | KIFAP3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | FLA3; KAP3; SMAP; KAP-1; KAP-3; Smg-GDS; dj190i16.1; KIFAP3                                                       |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 22920                                                                                                             |
| GeneID           | 22920                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human KIFAP3 (NP_055785.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, IP, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q92845                                                                                                            |
| Protein Weight   | 91kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |