

Anti-Kifc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bw19

Item Number: CM0029686



Introduction

KIFC1 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Targets human KIFC1 (HSET), a minus end-directed motor protein essential for spindle formation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	KIFC1 Rabbit pAb
Remarks/Alias	HSET; KNSL2; KIFC1
Species	Human
GeneID	3833
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BW19
Protein Weight	74kDa
Shipping	Ice bag