

Anti-Kifc1 Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - Q9Bw19

Item Number: CM0003410



Introduction

CamelBio supplies anti-KIFC1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse KIFC1. Essential for cell division and motor protein research.

[Learn More](#)

Parameter	Value
Product Name	KIFC1 Rabbit mAb
Remarks/Alias	HSET; KNSL2; KIFC1
Species	Human
Gene ID (Human)	3833
Gene ID	3833
Immunogen	Synthetic peptide within amino acids 1-100 of human KIFC1 (Q9BW19)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BW19
Protein Weight	74kDa
Shipping	Ice bag