

Anti-Kifc3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bvg8

Item Number: CM0030201



Introduction

Rabbit polyclonal antibody against human KIFC3, validated for WB and ELISA. Reacts with human and rat KIFC3, a minus-end microtubule motor protein involved in apically targeted transport and zonula adherens maintenance.

[Learn More](#)

Parameter	Value
Product Name	KIFC3 Rabbit pAb
Remarks/Alias	KIFC3
Species	Human
GeneID (Human)	3801
GeneID	3801
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 140-390 of human KIFC3 (NP_005541.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BVG8
Protein Weight	93kDa
Shipping	Ice bag