

Anti-Kinesin Light Chain 1 (Klc1) Polyclonal Antibody For Wb, Ihc-P, Elisa - Q07866

Item Number: CM0031062



Introduction

KLC1 Rabbit pAb detects human, mouse, rat KLC1 (UniProt Q07866) validated for WB, IHC-P, ELISA. A kinesin motor protein involved in organelle transport, ideal for neuroscience and intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	Kinesin light chain 1 (KLC1) Rabbit pAb
Remarks/Alias	KLC; KNS2; KNS2A; Kinesin light chain 1 (KLC1)
Species	Human
GeneID (human)	3831
GeneID	3831
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 261-560 of human Kinesin light chain 1 (KLC1) (NP_005543.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07866
Protein Weight	65kDa
Shipping	Ice bag