

Anti-Kinesin Light Chain 1 (Klc1) Rabbit Monoclonal Antibody For Wb, Elisa - Q07866

Item Number: CM0002446



Introduction

Rabbit monoclonal antibody targeting Kinesin light chain 1 (KLC1) for Western blot and ELISA. Cross-reacts with human and rat. Ideal for organelle transport and microtubule motor protein studies.

[Learn More](#)

Parameter	Value
Product Name	Kinesin light chain 1 (KLC1) Rabbit mAb
Remarks/Alias	KLC; KNS2; KNS2A; Kinesin light chain 1 (KLC1)
Species	Human
GeneID (human)	3831
GeneID	3831
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 450-550 of human Kinesin light chain 1 (KLC1) (Q07866)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q07866
Protein Weight	65kDa
Shipping	Ice bag