


Anti-Kininogen 1 (Kng1) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P01042

Item Number: CM0021422



Introduction

Rabbit polyclonal antibody against human Kininogen 1 (KNG1), validated for WB, IF/ICC, and ELISA applications; detects human, mouse, and rat KNG1. Suitable for coagulation and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	Kininogen 1 (KNG1) Rabbit pAb
Remarks/Alias	BK; HK; BDK; KNG; HAE6; HMWK; Kininogen 1 (KNG1)
Species	Human
Gene ID (Human)	3827
Gene ID	3827
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 148-427 of human Kininogen 1 (KNG1) (NP_000884.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01042
Protein Weight	72kDa
Shipping	Ice bag