


Anti-Kininogen 1 (Kng1) Rabbit Polyclonal Antibody For Wb - P01042

Item Number: CM0013444



Introduction

Polyclonal rabbit antibody targeting human Kininogen 1 (KNG1), validated for WB and ELISA. Detects human, mouse, rat kininogen; 72 kDa. Ideal for coagulation, bradykinin, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	Kininogen 1 (KNG1) Rabbit pAb
Remarks/Alias	BK; HK; BDK; KNG; HAE6; HMWK; Kininogen 1 (KNG1)
Species	Human
GeneID (human)	3827
GeneID	3827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 148-427 of human Kininogen 1 (NP_000884.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01042
Protein Weight	72kDa
Shipping	Ice bag