


Anti-Kir2DL2/Kir2DL3 Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P43628

Item Number: CM0012737



Introduction

Recombinant rabbit monoclonal antibody against human KIR2DL2/KIR2DL3, suitable for WB, IF/ICC, FC, and ELISA. Recognizes a natural killer cell inhibitory receptor for HLA-C ligands, expressed on NK cells.

[Learn More](#)

Parameter	Value
Product Name	KIR2DL2/KIR2DL3 Rabbit mAb
Remarks/Alias	KIR2DL3; CD158B2; CD158b; GL183; KIR-023GB; KIR-K7b; KIR-K7c; KIR2D55; KIRCL23; NKAT; NKAT2; NKAT2A; NKAT2B; p58; killer cell immunoglobulin-like receptor 2DL3; KIR2DL2/KIR2DL3
Species	Human
Gene ID (Human)	3803/3804
Gene ID	3803/3804
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human KIR2DL3(NP_056952.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P43628
Protein Weight	27kDa/37kDa/38kDa
Shipping	Ice bag