

Anti-Kir2DL4 Monoclonal Antibody For Wb, Fc, Elisa - Q99706

Item Number: CM0007421



Introduction

Rabbit monoclonal antibody against human KIR2DL4, validated for Western blot, flow cytometry, and ELISA. This receptor binds HLA-G and regulates immune responses at the maternal-fetal interface.

[Learn More](#)

Parameter	Value
Product Name	KIR2DL4 Rabbit mAb
Remarks/Alias	G9P; CD158D; KIR103; KIR-2DL4; KIR103AS; KIR-103AS; KIR2DL4
Species	Human
GeneID (Human)	3805
GeneID	3805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-242 of human KIR2DL4 (NP_002246.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	Q99706
Protein Weight	41kDa
Shipping	Ice bag