


Anti-Kir3DL1 Polyclonal Antibody For Wb, Elisa - P43629

Item Number: CM0023405



Introduction

Anti-KIR3DL1 polyclonal antibody validated for WB and ELISA. Targets human Killer cell immunoglobulin-like receptor 3DL1, cross-reacts with mouse, MW 49 kDa. Ideal for NK cell and immunoreceptor studies.

[Learn More](#)

Parameter	Value
Product Name	KIR3DL1 Rabbit pAb
Remarks/Alias	KIR; NKB1; NKAT3; NKB1B; NKAT-3; CD158E1; KIR2DL5B; KIR3DL1/S1; KIR3DL1
Species	Human
Gene ID (Human)	3811
Gene ID	3811
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 22-222 of human KIR3DL1 (NP_037421.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P43629
Protein Weight	49kDa
Shipping	Ice bag